

Centre for Advanced Electron Spin Resonance

"a state-of-the-art facility for use by
biochemists, chemists and physicists"


Dr Will Myers, 2 October, 2025, ICL Lecture Theatre

Introductory Lecture

- User Applications
- Instrumentation
- ESR Personnel
- How to use CAESR
- Magnet Hazards, Chemical,
- Microwave, and Laser Safety
- Sample Preparation
- Data Storage & Processing
- User Resources



Clarendon




ICL building

User Applications

- Diamagnetic \rightarrow Paramagnetic
 - 20 Å to 120 Å
 - Synthesize & Characterize labels:
 - Nitroxides, Trityl, Gd(III);
 - Redox/ E-chem to open shell states
- Intrinsic Paramagnets
 - Paramagnetic Catalyst Intermediates
 - Metallo-enzyme Mechanisms
 - Single Molecule Magnetism
- Transient Paramagnets
 - Photo-Activated Transient States
 - 1 ns resolution
 - Electrical &/or Optical Det

Pulsed EDMR






CW-EPR Characterisation

$$\gamma(e^-) = -658 * \gamma(^1H)$$

DEPARTMENT OF
CHEMISTRY





UNIVERSITY OF
OXFORD



$$h\nu = g_{eff} \mu_B B_0$$

$$g_{eff} = \sqrt{g_{||}^2 \cos^2 \theta + g_{\perp}^2 \sin^2 \theta}$$

neat solid, room temperature




JMG DFT predicts $g_{||} = 2.000$ $g_{\perp} = 2.023$

by experiment, $g_{||} = 1.993$ $g_{\perp} = 2.043$

Magnetic dilution, ← Not just e⁻ hyperfine resolution

What about ⁷³Ge? I = 9/2

JMG DFT predicts 4x ax. hfc 3 MHz
8x eq. hfc -20 MHz




$$\text{Pr}(k) = n! f^k (1-f)^{n-k} / k!(n-k)!$$


n total number of sites

k number of occupied sites

f ⁷³Ge abundance, 7.76%






3 mM in
pyridine,
100K



Hyperfine Resolution, Beyond CW-EPR

NiFe Hydrogenase



Dark-adapted sample,
mixture of NiC & NiL

NiL

30 min. of laser pulses
1 mJ pulse energy
550 nm wavelength
20 Hz repetition rate




NiC

90 min. annealing
at 200K

NiC simulation

$A(^1\text{H}) = [18.4 \ -10.8 \ -18] \text{ /MHz}$
Lit. ENDOR of *R. eutropha*,
Brecht, *et al.*, JACS 2003, 125, 13075.

HYSCORE measured at 2.5 K



ESR Instrumentation

5 research instruments

<u>name</u>	<u>GHz</u>	<u>band</u>	<u>method</u>	<u>location</u>
EMX _{MICRO}	9.1 - 9.9	X	CW	ICL F19
E580	9.1 - 9.9	X	CW / Pulsed	ICL F11
	33 - 35	Q	CW / Pulsed	
E680	9.1 - 9.9	X	CW / Pulsed	ICL F12
	92 - 94	W	CW / Pulsed	
E380	9.1 - 9.9	X	CW / Pulsed	Clarendon 020
Krymov	130	D	CW / Pulsed	Clarendon 020

EMX_{MICRO}


ICL/F19

- CW-EPR only
- Excellent SNR
- 2.5 – 300 K with ESR-900
- 100 - 450 K with N₂ heater
- 77 K with N_{2(l)} finger dewar
- Automatic goniometer
- Room interlocked for Class 4/3B lasers
- User scheduling



E680


ICL/F12



- 9.1-9.9 GHz, X-band & 94 GHz, W-band
- 2.5 – 300 K with CF-935 & ESR-900
- OPO laser: 213-1700 nm, 6-100 mJ, 7 ns pulse length 20 Hz
- API & ProDel automation
- Arbitrary Waveform Generator (AWG)
- Room interlocked for Class 4/3B lasers

E580

ICL/F11



- X-band and Q-band
- Arbitrary Waveform Generator (AWG)
- 140 W at Q-band
- ProDel & API - programmable
- 2.5 – 300 K with CF-935 & ESR-900 cryostats
- OPO laser 355, 410-2200 nm; 4-13 mJ; 7 ns length, 20 Hz
- Room interlocked for Class 4/3B lasers

CAESR Personnel



Professor Christiane Timmel
<http://research.chem.ox.ac.uk/christiane-timmel.aspx>

Professor Arzhang Ardavan
<https://www2.physics.ox.ac.uk/contacts/people/ardavan>



Dr Will Myers



Office 72651
ICL F20a

Lab 72665
ICL F12

Fax 75410

william.myers@chem.ox.ac.uk

<http://research.chem.ox.ac.uk/william-myers.aspx>

1. Run CAESR & Services in Chemistry
4. Consultative role to students and staff
5. Collaborate with Department and External groups
6. Assist in joint publications
7. Pursue, publish, and present peer-reviewed research
9. Involved in the submission of grant applications
14. Students and staff training
15. Responsible for facility safety regulations
- Assisting w/ supervision of Part II and D. Phil. student projects
- Booking Meetings
- Teaching: Lecturing, small-group teaching, tutoring

How to Use the CAESR SRF

1. **Independent:** Previous ESR (EPR) or support of experienced group member, ask for any help.
2. **Semi-independent:** Complete training, then independent measurements, scheduling, and data analysis. User asks for support as required.
3. **With CAESR staff:** Complete training, then semi-independent measurements and scheduling. Extended support with project development, data analysis and publications.
4. **Service:** “One-off” projects are rare. Routine or repetitive measurements are generally not possible, (see #2).
5. **Collaboration:** Collaborate with the Scientific Applications Manager, or a D Phil from one of the co-Director groups might be available to use your project for their degree course.

Typical User Sequence for ESR

- Discussion with Scientific Applications Manager
 - Research plan, safety, and literature
 - Preliminary sample list
 - Purchase of clear-fused Quartz EPR tubes
- Booking CW-EPR time, COSHH
- Sample Preparations in F13
- CW-EPR in F19, with any required training
- Monthly booking meeting for Pulsed EPR
- Pulsed EPR and / or W-band CW-EPR

Magnetic Field Safety in ESR

Strong Magnetic Fields


Hazards to:

- heart pacemakers
- magnetic back or ID cards
- watches (non-LCD)

Stray fields in Corridors

- outside F12 at a max field

Electromagnet field is mostly
within its yoke




Cryogenic & Chemical Safety in ESR

- ❖ Wet lab work in F13 → departmental guidelines
- ❖ All samples (hazardous or not), need COSHH form on file in F20a.
- ❖ Samples occasionally break -- act quickly to save before complete thawing.
- ❖ Safety courses are available for how to handle liquid nitrogen and liquid helium.

LASER Safety in ESR

- ☐ All CAESR labs are inter-locked for use of Class 4 and 3B lasers.
- ☐ Pulsed OPO lasers for 230 - 2200 nm, 4 - 100 mJ, 20 Hz available for Transient ESR. CW lasers at 355 and 532 nm are available.
- ☐ New laser users must attend Chemistry Department safety lecture
- ☐ Users who wish to use their own lasers in CAESR must follow Departmental laser installation procedures.
- ☐ Laser Safety forms for each lab must be signed.

Microwave Safety in ESR



$\lambda \sim 3.2 \text{ cm}$

- Place microwave sources in Standby when changing resonators and samples.
- DON'T stare into the end of an open waveguide with microwaves propagating outward.
- X-band (9.4 GHz) exposure may result in cataracts.
 - ✓ approx. 60 min of CW at 200 mW / cm²

CAESR Wet Lab ICL/F13



- Drying oven, balance, sonicator, vortex mixer
- Fume Cupboard and Schlenk/Vacuum Line
- Flame-sealing tubes after de-gassing via FREEZE-PUMP-THAW
- Flame-sealing J. Young tube bottoms
- Mbraun LABstar ; Cary60 UV-vis

Sample tubes and volumes

Tube Dimensions			Micro-wave Band	Sample Vol. (μL)	<u>Comment</u>
O.D. (mm)	I.D. (mm)	L (mm)			
0.84	0.7	100	W	3	General, for powders and frozen solutions
1.6	1.2	100	Q	10	Dielectric ENDOR resonator
3	2.6	150	X	40	For split ring resonator
3.80	2.79	150	X	80	For ENDOR pulsed resonator, high precision wall thickness
4	3	250	X	250	For all CW resonators and MD5, low-dielectric solutions/solids only

Tubes of are for sale in ICL/F20a, and on the R12 Oracle purchasing system.

Sample Concentrations


Transition metals 1 - 5 mM

Organic radicals 20 - 200 mM

Solids 1×10^{15} to 1×10^{18} spins

Sample concentration vs. *Spin* concentration

Single Crystals – 1×10^9 – auto goniometer



Equal concentrations, Equal noise,
Different ESR signal width

Data Storage and Transfer

Storage

On the Spectrometers: Data is kept in
.../xeprfiles/data/YourPI_group/You/DateName/...

Back-Up: \\chem.ox.ac.uk\SRF\ESR\Spectrometer_DataBackup\...

→ But no guarantees, you are responsible for your data.

Transfer

File Transfer Protocol (FTP) is available on E580, E680, and EMXmicro

See: \\chem.ox.ac.uk\SRF\ESR\
NewUser_IntrosAndFAQ\FTP_FileTransfer_conversion.pdf

Data Processing Software

Windows, Mac and Linux:

Matlab - dept. - Data plotting, integration, spin quantitation, simulations

Python, C, C++, Java, Fortran, MS Excel, etc. -free -

EasySpin – free, w/ Matlab - www.easyspin.org – simulation software

Spinach – free, w/ Matlab - www.spindynamics.org - simulation software

SpinDynamica – free, w/ Mathematica - www.spindynamica.org -
simulation software

+ many others

Windows:

SpinCount - \$ - Prof. M. Hendrich, Carnegie Mellon Univ. (not in CAESR)

Linux:

XSophe - \$\$\$\$ - Bruker BioSpin (not in CAESR)

User Resources

3rd year EPR practical

See: [\\chem.ox.ac.uk\\SRF\\ESR](http://chem.ox.ac.uk/SRF/ESR)

Video lecture links

Lecture course slides

Spectrometer Manuals

Stop by -or- questions by e-mail

- Main Reference Books -- RSL
 - Pulse EPR methods
 - Selected Systems
 - Instrumentation
- Conferences
- Organizations
- Databases and Software
- Spectrometer Makers
- List of 169 Reference Books

<https://caesr.web.ox.ac.uk/home>

Graduate Training Sessions

Morning ~ 10-12 pm

<https://forms.office.com/e/zGrkuguJFQ>

Explanation of webpage and online resources

Tour of ESR Facility

How to Set up cryogenics

How to use a Spectrometer

Send an interest e-mail to
william.myers@chem.ox.ac.uk

Instrumentation description

Demonstration: CW-EPR & Pulsed EPR of LiPc

Pulsed ESR of BDPA: Electron-Nuclear DOuble Resonance

Afternoon ~ 2:30 to 5:00 pm

Demonstration: Transient ESR *with* LASER excitation

Double Electron-Electron Resonance

Bruker Pulsed EPR Course – group use of the spectrometer